

# A-Flex™ Plus Protective Headgear



Orthomerica introduces the A-Flex Plus Protective Headgear, an adjustable alternative to other hard protective devices. The A-Flex Protective Headgear is fabricated from a flexible plastic that readily conforms to varying head shapes making it easy-to-fit upon demand. The protection is ideal for low impact forces that are distributed across the entire surface without cracking or penetrating. The headgear is equipped with chin straps that are recommended for active and noncompliant patients. The posterior strap “locks” the device below the occiput and mastoids and accommodates size variations. The lateral straps are used for easy donning and doffing. Lateral straps and side-slits are also optional. The A-Flex features ventilation holes to improve air circulation and cooling. It can be easily trimmed with a pair of heavy duty scissors without the need for complex equipment and tools. It can be stocked in various sizes.

## Clinical Indications

- Post-operative craniectomy or craniotomy
- Cranial vault reconstruction or other cranial surgery
- Skull anomalies
- Seizure disorders
- Self-injurious behaviors related to head banging, hair pulling, etc.
- Instability secondary to poor balance and/or disequilibrium reactions

## Contraindications

- Individuals who are at high risk for facial injuries (Custom products including face guards and other forms of cranial protection are available from Orthomerica)
- Individuals requiring high impact protection

Size	Blue	White	Circumference
Pediatric SM	<b>783.805</b>	<b>783.905</b>	44-47 cm
Pediatric MD	<b>783.815</b>	<b>783.915</b>	46-49 cm
Pediatric LG	<b>783.825</b>	<b>783.925</b>	48-51 cm
Pediatric XL	<b>783.835</b>	<b>783.935</b>	50-53 cm
Adult XS	<b>783.845</b>	<b>783.945</b>	52-55 cm
Adult SM	<b>783.855</b>	<b>783.955</b>	54-57 cm
Adult MD	<b>783.865</b>	<b>783.965</b>	57-60 cm
Adult LG	<b>783.875</b>	<b>783.975</b>	60-63 cm
Adult XL	<b>783.885</b>	<b>783.985</b>	63-66 cm
Adult 2X	<b>783.895</b>	<b>783.995</b>	66-69 cm

Chin Strap Kit—Blue	<b>783.15</b>
Chin Strap Kit—White	<b>783.13</b>
Chin Strap Kit—Transfer Pattern	<b>783.16</b>
Replacement Side Strap Kit—Blue	<b>783.17</b>
Replacement Side Strap Kit—White	<b>783.20</b>
Replacement Side Strap Kit—Custom	<b>783.19</b>

A-Flex Plus Protective Headgear may be ordered without the Chin Strap by excluding the “S” in the Model Number.  
Chin Straps are one-size-fits-all.

**ORTHOMERICA®**

(800) 446-6770 | [www.orthomerica.com](http://www.orthomerica.com)

6333 North Orange Blossom Trail • Orlando, FL 32810 • USA



Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands