

ORTHOMERICA® OWLS® WHO™ & shark-o™

6333 North Orange Blossom Trail, Orlando FL 32810 • www.orthomerica.com • phone (877) 737-8444 • fax (877) 737-8445

Patient last name/ID	Age	Weight <small>(if weight is 220 lbs. [100 kg] or more, ankle must be reinforced at added cost)</small>	Activity level	Edema present?	Shipping Provider UPS, FedEx, Other Service Ground, 2 day air, 3 day air, Overnight
Patient first name	Height	Low Medium High <small>If high, use one weight range higher.</small>		None Light Medium Heavy	
Male, Female	Diagnosis				
Facility	Date	PO#			Need by
Practitioner	Ship-to address			Billing address	
Email	City	State	City	State	

Phone	Fax	Zip	Country	Zip	Country
Left, Right <small>(If bilateral use two forms For bilateral CLAW, use one weight range higher)</small>					
Cast, Scan <small>(Tracings are recommended for cast orders)</small>	shark-o <small>(charcot orthosis)</small>	WHO Heel <small>(cutout relief)</small>	WHO Forefoot <small>(cutout relief)</small>	WHO Midfoot <small>(full foot bed)</small>	WHO Chopart <small>(partial foot)</small>
					CLAW™ Carbon Laminate Alignment Walker <small>(available on all WHO and shark-o models, added cost, patten bottom, common HCPCS code L-2370)</small>

For Chopart level, use one weight range higher.
Over "toe" strap required (to keep brace from flaring).

Ankle flexion
 As casted (default)
 90°
 Other

Forefoot
 As casted (default)
 Neutral
 Other

Hindfoot
 As casted (default)
 Subtalar neutral
 Other

Liner
 Black (default)
 White

Copolymer
 1/4" (default)
 3/16"


Copolymer color
 Black (default)
 Natural

Pattern transfer
(added cost, over natural plastic only, available patterns on www.orthomerica.com)

Carbon ankle reinforcements
 Select to add reinforcements
(added cost, required if pt. weight 220 lbs. or above or if CLAW is ordered.)

Measurements

If anticipated CLAW alignment more than 20 degrees, how much in or out?



- Circumferences**
 A Calf
 B Forefoot
 C Ankle
- Lengths**
 F Brace Height
 G Foot Length
- Diameter**
 D Forefoot ML
 E Ankle ML

Mark ulcer/amputation sites
(draw in Adobe Acrobat: select Comment, Drawing Markups, Draw Free Form)



Foot tracing recommended

Evenup contralateral height compensation
(added cost, slip-on, over the shoe sandal, with two each 1/2" build-up filler insoles)



Enter shoe size to order

Foot weight bearing during measurement? Yes, No

Remarks