

Orlando[®] KAFO/HKAFO

LOWER EXTREMITY FRACTURE BRACING SYSTEM

For over 20 years, Orthomerica Products, Inc. has been offering effective non-operative and post-operative treatment solutions for lower extremity fractures and range-of-motion restriction.

Orlando KAFO/HKAFOs are modular, lightweight and can be ordered and stocked as components, making them immediately available, or custom assembled to your measurements. Orthotists choose from a wide selection of femoral and tibial components, hip joints, knee joints and shoe inserts to meet individual patient needs.

Anatomically designed pelvic components offer a variety of shapes and sizes to accommodate even the most difficult to fit patient. Unique, patented hip joints control flexion/extension and abduction/adduction. These joints control range of motion of the affected limb and provide adjustability throughout the rehabilitation process. Thigh and tibial components promote further stability through circumferential compression of soft tissue. Standard or special order knee joints control sagittal plane motion. The shoe insert provides the final link in the mechanical chain of support by allowing limited or full range of ankle motion as indicated.


INDICATIONS

- Post-op hip revision patients
- Post-op lower extremity fractures
- Patients needing stability following the removal of external fixation devices
- Non-operative fractures
- Fracture management following internal fixation
- Orthotic rotational control for patients with chronic hip dislocation

FEATURES

- HKAFO limits hip rotation and provides adjustable abduction/adduction and flexion/extension control
- Choice of several tibial components, hip joints, knee joints and shoe components
- Lined or unlined femoral/tibial components
- Custom-to-measurements phoned, faxed or e-mailed for same day shipment
- Prefabricated components can be assembled on-site

Available in individual prefabricated components, pre-assembled, custom-to-measurements (all shipped same day) or custom-to-cast.



ORTHOMERICA® CTM ORLANDO KAFO/HKAFO ORDER FORM

6333 North Orange Blossom Trail, Orlando FL 32810 • www.orthomerica.com • phone 800-446-6770 • fax 800-638-9259

PATIENT INFO (PH)

LAST NAME/ID _____ FIRST NAME _____

AGE _____ HEIGHT _____ WEIGHT _____ SEX _____

DIAGNOSIS _____

PO# _____

REMARKS: _____

SHIPPING INFORMATION

PRACTITIONER _____ PHONE _____ FAX _____

FACILITY _____

SHIP TO ADDRESS _____


CITY _____ STATE _____ ZIP _____

Ground 2 Day Air Overnight Other: _____

Need by: _____

BILL TO ADDRESS _____

CITY _____ STATE _____ ZIP _____



- A = Waist Circumference
- B = Hip Circumference
- C = Proximal Thigh Circ.
- D = Distal Thigh Circ.
- E = Largest Calf Circ.
- F = (Just above) Ankle Circ.
- G = Length, Perineum to Floor
 .1 Anatomical
 .2 Actual Medial Brace
- H = Caliper Meas., Knee Center (M/L)
- I = Length, Knee Center to Floor
- J = Length, Trochanter to Knee Center
- K = Waist to Trochanter

STYLE

HKAFO KAFO
(KAFO only requires dimensions C-I)

Left Right

PELVIC SECTION

Newport® 3 Newport® 4
 Newport® MC (CTM 2 strap upgrade charge)
 Newport® Jr.

HIP JOINTS

Virtual® 5
 Virtual® (adj. ab/adduction, flexion/extension)

Following WILL NOT work with Newport 4:
 Adj. ROM (adj., flexion/extension, 15° abduction)
 Drop Lock*

THIGH COMPONENTS

Ischial NSNA

KNEE JOINTS

(Special order knee joints available—please call for more info)

Free Motion (1/8" 3/16"*)
 Drop Lock (1/8" 3/16"*)
 Adj. ROM* (1/8" 3/16")

TIBIAL COMPONENTS

Orlando Tibial
 Orlando Lite Tibial section
(with bilateral uprights)
 Orlando TPFx

SHOE INSERTS

Free Motion
 Adjustable Single Axis
(Extended footplate, locks at 90° or has free ROM)

Solid Ankle
 PRAFO*

LINER*



Aliplast (1/8" 3/16" 1/4")
 CoolFoam®

FRACTURE SOCKS (adult)

Yes, Quantity _____
 No

*** Incurs upgrade charge**

Hip Joints **Knee Joints**



Emergo Europe
 Molenstraat 15, 2513 BH
 The Hague, The Netherlands
 Phone +31-70-345-8570
 Fax +31-70-346-7299



6333 N Orange Blossom Trail
 Orlando FL 32810 USA
 800-446-6770
 orthomerica.com