

New!

California[®] High Profile TLSO™



The California High Profile TLSO is available in two designs (460/464) and requires customization made by a qualified clinician (orthotist) to effectively fit each patient. The High Profile TLSO provides stability and triplanar control where indicated to reduce spinal muscle atrophy and deconditioning.

Features

- Patented compound closure system provides patient-regulated compression to supplement abdominal musculature
- Coolfoam™ liner provides a comfortable patient interface
- Thermoplastic posterior component extends to T9 and provides total contact over the paraspinal muscles
- Anterior panel includes a pocket with an internal thermoplastic shell to provide abdominal support and control motion
- Control panels provide additional lateral strength in 464 model

Indications

- Post-operative laparoscopic disc replacement
- Post-operative IDET procedures (intradiscal electrothermal annuloplasty)
- Post-operative lumbar laminectomy
- Post-operative spinal fusion
- Chronic low back pain
- Mechanical back pain with activities of daily living
- When spinal control of extension and flexion is desired/required
- Compression fractures

Size	Fits Waists	Standard ATE		Extended ATE	
XS	25"-30"	3239.460.01	3239.464.01	3239.460.02	3239.464.02
SM	30"-35"	3240.460.01	3240.464.01	3240.460.02	3240.464.02
MD	35"-40"	3241.460.01	3241.464.01	3241.460.02	3241.464.02
LG	40"-45"	3242.460.01	3242.464.01	3242.460.02	3242.464.02
XL	45"-50"	3243.460.01	3243.464.01	3243.460.02	3243.464.02
2XL	50"-55"	3244.460.01	3244.464.01	3244.460.02	3244.464.02
3XL	55"-60"	3245.460.01	3245.464.01	3245.460.02	3245.464.02



**6333 North Orange Blossom Trail
Orlando FL 32810**

(800) 446-6770 | www.orthomerica.com

