


CUSTOM -TO- MEASUREMENTS

NEWPORT[®] HIP SYSTEM



MEASUREMENT LEGEND


A = Circumference at Xyphoid Process
 B = Circumference at Waist
 C = Circumference at Hip
 D = Waist to Trochanter
 E = Proximal Thigh Circumference
 F = Distal Thigh Circumference
 G = Hip to Knee Center
 H = Knee Center to Floor

- Black CoolFoam™ Liner *
- Hip Interface Sock *
_____ Qty. (1 per pkg.)
- Kodel Liner *

PELVIC BAND STYLE :

- Maximum Control (2 straps/CTM only)
- Maximum Control (2 straps, high contralateral, 1 side only)
- Newport³ (3" pelvic height only)
- Newport⁴
- Newport^{jr}

- ARC™ Kit ****
- ARC™ Kit with PRAFO ****



Shoe Size _____

HIP JOINT OPTIONS :

Standard Virtual® Hip joint will be supplied unless otherwise specified

Virtual® Hip Joint (adjustable abduction/adduction, flexion/extension)



- Lerman style *** (adjustable flexion/extension, 15° abduction)
- * Drop Lock *** (15° abduction)
- * Chronic Joint *** (adjustable abduction, flexion/extension)
- NP jr. Abduction Bar

NEWPORT® 4 HIP JOINT OPTIONS :
(Available on Newport 4 pelvic only)



- Virtual® 5 For anterior or posterior dislocations and fractures

* Additional Upgrade Charge
 ** Additional Upgrade Charge for:
 circumferences > 45"
 circumferences > 60"
 *** Not Available with NP jr.

ORDER INFORMATION Account #: _____

Date: _____ P.O.#: _____

Facility to be Billed: _____

Ship to Address: _____

Phone: _____ Fax: _____

Date Required: _____

Ship Via: _____ On (Date) _____

Contact: _____

PATIENT INFORMATION (PHI)

Date: _____

Patient ID/Name: _____

Age: _____ Sex: _____ Height: _____ Weight: _____

Diagnosis: _____

**Fax This Form to 800.638.9259
 or Call 800.446.6770**

U.S. Patent No. 5,421,810

SL00100 REV A
 © 2008 ORTHOMERICA PRODUCTS, INC. ALL RIGHTS RESERVED.

E-mail: custserv@orthomerica.com
<http://www.orthomerica.com>

ORTHOMERICA[®]
 THE GLOBAL ORTHOTIC SOLUTIONSM