

**NEW!**

READY-TO-FIT PREFABRICATED LSO/TLSO

# STOCK SPINAL ORTHOSES

YOUR SOLUTION for  
**5pm Friday  
& Weekend**  
patients.

Stock Spinal LSO/TLSO fits  
75% of Adult Patient  
Population



**Bivalve LSO with Corset Front and PSR (Posterior Spinal Relief) Back**



**Bivalve TLSO with ATE (Anterior Thoracic Extension)**



**California Soft Spinal Front TLSO with ATE**

## Features and Benefits

- Keep inventory in your office for same day fitting
- Initial stocking order of six or more braces receives discount
- Have inventory on hand for late night/weekend fittings
- Can reduce the need for 24-hour hospital stays
- Save cost of overnight shipping
- Braces are finished and ready to fit
- Twelve sizes—six male, six female (XS-XXL) in neutral or P.A. fronts in 0° or 15° lordosis
- Available as an LSO, without ATE
- Extended ATE available

## Standard Materials

- $\frac{5}{32}$ " Mod Natural Polyeth
- $\frac{3}{16}$ " Aliplast™, Volara®
- Adjustable 1½" Dacron Straps with Chicago Screws
- Straps Anterior
- Extended Chafes

**ORTHOMERICA®**

# STOCK SPINAL CORSET FRONT LSO

Order Number:

## Customer Information

Date Company PO# Email

Contact Phone Fax

Address City State Zip

Ship Via/Shipping Instructions Date Required

### Women's Corset Front with 15° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	2212.01	28"-31"	29"-32"	25"-28"	34"-37"
	Small	2212.02	30"-34"	31"-35"	27"-31"	36"-39"
	Medium	2212.03	32"-36"	34"-38"	30"-34"	38"-42"
	Large	2212.04	37"-41"	37"-41"	35"-39"	41"-45"
	X-Large	2212.05	40"-44"	41"-45"	40"-44"	42"-48"
	2X-Large	2212.06	44"-48"	46"-50"	45"-49"	50"-54"

### Women's Corset Front with 0° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	2212.13	28"-31"	29"-32"	26"-29"	34"-38"
	Small	2212.14	30"-34"	31"-35"	28"-32"	36"-39"
	Medium	2212.15	32"-36"	34"-38"	31"-35"	38"-42"
	Large	2212.16	37"-41"	37"-41"	36"-40"	41"-45"
	X-Large	2212.17	40"-44"	41"-45"	41"-45"	43"-49"
	2X-Large	2212.18	44"-48"	46"-50"	46"-50"	50"-54"

### Men's Corset Front with 15° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	2212.07	34"-38"	35"-39"	30"-34"	36"-40"
	Small	2212.08	36"-40"	36"-40"	33"-37"	38"-42"
	Medium	2212.09	39"-43"	38"-42"	36"-40"	42"-46"
	Large	2212.10	42"-46"	42"-44"	40"-44"	44"-48"
	X-Large	2212.11	45"-49"	46"-50"	45"-49"	47"-51"
	2X-Large	2212.12	50"-54"	52"-56"	50"-54"	53"-57"

### Men's Corset Front with 0° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	2212.19	34"-38"	35"-39"	31"-35"	36"-40"
	Small	2212.20	36"-40"	36"-40"	34"-38"	38"-42"
	Medium	2212.21	39"-43"	38"-42"	37"-41"	42"-46"
	Large	2212.22	42"-46"	42"-44"	41"-45"	44"-48"
	X-Large	2212.23	45"-49"	46"-50"	46"-50"	47"-51"
	2X-Large	2212.24	50"-54"	52"-56"	51"-55"	53"-57"

## Customer Information

Date Company PO# Email

Contact Phone Fax

Address City State Zip

Ship Via/Shipping Instructions Date Required

### Women's Neutral Abdomen with 15° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2207.01</b>	28"-31"	29"-32"	25"-28"	34"-37"
	Small	<b>2207.02</b>	30"-34"	31"-35"	27"-31"	36"-39"
	Medium	<b>2207.03</b>	32"-36"	34"-38"	30"-34"	38"-42"
	Large	<b>2207.04</b>	37"-41"	37"-41"	35"-39"	41"-45"
	X-Large	<b>2207.05</b>	40"-44"	41"-45"	40"-44"	42"-48"
	2X-Large	<b>2207.06</b>	44"-48"	46"-50"	45"-49"	50"-54"

### Women's Pendulous Abdomen with 15° Lordosis

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2207.01PA</b>	28"-32"	29"-33"	26"-30"	35"-39"
	<b>2207.02PA</b>	30"-34"	31"-35"	29"-33"	37"-41"
	<b>2207.03PA</b>	32"-36"	34"-38"	32"-36"	38"-42"
	<b>2207.04PA</b>	36"-40"	37"-41"	37"-41"	41"-45"
	<b>2207.05PA</b>	41"-45"	42"-46"	43"-47"	45"-49"
	<b>2207.06PA</b>	43"-47"	46"-50"	45"-49"	50"-54"

### Women's Neutral Abdomen with 0° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2207.13</b>	28"-31"	29"-32"	26"-29"	34"-38"
	Small	<b>2207.14</b>	30"-34"	31"-35"	28"-32"	36"-39"
	Medium	<b>2207.15</b>	32"-36"	34"-38"	31"-35"	38"-42"
	Large	<b>2207.16</b>	37"-41"	37"-41"	36"-40"	41"-45"
	X-Large	<b>2207.17</b>	40"-44"	41"-45"	41"-45"	43"-49"
	2X-Large	<b>2207.18</b>	44"-48"	46"-50"	46"-50"	50"-54"

### Women's Pendulous Abdomen with 0° Lordosis

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2207.13PA</b>	28"-33"	30"-34"	30"-34"	35"-39"
	<b>2207.14PA</b>	31"-35"	32"-36"	32"-36"	36"-40"
	<b>2207.15PA</b>	33"-37"	35"-39"	35"-39"	39"-43"
	<b>2207.16PA</b>	37"-41"	38"-42"	38"-42"	41"-45"
	<b>2207.17PA</b>	42"-45"	43"-47"	43"-47"	45"-49"
	<b>2207.18PA</b>	44"-48"	48"-52"	49"-53"	51"-54"

### Men's Neutral Abdomen with 15° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2207.07</b>	34"-38"	35"-39"	30"-34"	36"-40"
	Small	<b>2207.08</b>	36"-40"	36"-40"	33"-37"	38"-42"
	Medium	<b>2207.09</b>	39"-43"	38"-42"	36"-40"	42"-46"
	Large	<b>2207.10</b>	42"-46"	42"-44"	40"-44"	44"-48"
	X-Large	<b>2207.11</b>	45"-49"	46"-50"	45"-49"	47"-51"
	2X-Large	<b>2207.12</b>	50"-54"	52"-56"	50"-54"	53"-57"

### Men's Pendulous Abdomen with 15° Lordosis

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2207.07PA</b>	34"-38"	34"-38"	32"-36"	35"-39"
	<b>2207.08PA</b>	36"-40"	36"-40"	35"-39"	38"-42"
	<b>2207.09PA</b>	39"-43"	38"-42"	38"-42"	41"-45"
	<b>2207.10PA</b>	42"-46"	42"-46"	42"-46"	44"-48"
	<b>2207.11PA</b>	45"-49"	47"-51"	47"-51"	48"-52"
	<b>2207.12PA</b>	51"-55"	51"-55"	51"-55"	53"-57"

### Men's Neutral Abdomen with 0° Lordosis

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2207.19</b>	34"-38"	35"-39"	31"-35"	36"-40"
	Small	<b>2207.20</b>	36"-40"	36"-40"	34"-38"	38"-42"
	Medium	<b>2207.21</b>	39"-43"	38"-42"	37"-41"	42"-46"
	Large	<b>2207.22</b>	42"-46"	42"-44"	41"-45"	44"-48"
	X-Large	<b>2207.23</b>	45"-49"	46"-50"	46"-50"	47"-51"
	2X-Large	<b>2207.24</b>	50"-54"	52"-56"	51"-55"	53"-57"

### Men's Pendulous Abdomen with 0° Lordosis

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2207.19PA</b>	34"-38"	35"-39"	34"-38"	36"-40"
	<b>2207.20PA</b>	36"-40"	38"-42"	38"-42"	39"-43"
	<b>2207.21PA</b>	37"-41"	40"-44"	41"-45"	42"-46"
	<b>2207.22PA</b>	40"-44"	44"-48"	44"-48"	44"-48"
	<b>2207.23PA</b>	45"-49"	48"-52"	49"-53"	47"-51"
	<b>2207.24PA</b>	50"-54"	53"-56"	54"-58"	53"-57"

# STOCK SPINAL CA SOFT SPINAL FRONT TLSO

Order Number:

**Customer Information**  No ATE (add an "L" to the part number)

Date Company PO# Email

Contact Phone Fax

Address City State Zip

Ship Via/Shipping Instructions Date Required

**Women's Neutral Abdomen with 15° Lordosis**

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2215.01</b>	28"-31"	29"-32"	25"-28"	34"-37"
	Small	<b>2215.02</b>	30"-34"	31"-35"	27"-31"	36"-39"
	Medium	<b>2215.03</b>	32"-36"	34"-38"	30"-34"	38"-42"
	Large	<b>2215.04</b>	37"-41"	37"-41"	35"-39"	41"-45"
	X-Large	<b>2215.05</b>	40"-44"	41"-45"	40"-44"	42"-48"
	2X-Large	<b>2215.06</b>	44"-48"	46"-50"	45"-49"	50"-54"

**Women's Neutral Abdomen with 0° Lordosis**

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2215.13</b>	28"-31"	29"-32"	26"-29"	34"-38"
	Small	<b>2215.14</b>	30"-34"	31"-35"	28"-32"	36"-39"
	Medium	<b>2215.15</b>	32"-36"	34"-38"	31"-35"	38"-42"
	Large	<b>2215.16</b>	37"-41"	37"-41"	36"-40"	41"-45"
	X-Large	<b>2215.17</b>	40"-44"	41"-45"	41"-45"	43"-49"
	2X-Large	<b>2215.18</b>	44"-48"	46"-50"	46"-50"	50"-54"

**Men's Neutral Abdomen with 15° Lordosis**

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2215.07</b>	34"-38"	35"-39"	30"-34"	36"-40"
	Small	<b>2215.08</b>	36"-40"	36"-40"	33"-37"	38"-42"
	Medium	<b>2215.09</b>	39"-43"	38"-42"	36"-40"	42"-46"
	Large	<b>2215.10</b>	42"-46"	42"-44"	40"-44"	44"-48"
	X-Large	<b>2215.11</b>	45"-49"	46"-50"	45"-49"	47"-51"
	2X-Large	<b>2215.12</b>	50"-54"	52"-56"	50"-54"	53"-57"

**Men's Neutral Abdomen with 0° Lordosis**

Quantity	Size	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	X-Small	<b>2215.19</b>	34"-38"	35"-39"	31"-35"	36"-40"
	Small	<b>2215.20</b>	36"-40"	36"-40"	34"-38"	38"-42"
	Medium	<b>2215.21</b>	39"-43"	38"-42"	37"-41"	42"-46"
	Large	<b>2215.22</b>	42"-46"	42"-44"	41"-45"	44"-48"
	X-Large	<b>2215.23</b>	45"-49"	46"-50"	46"-50"	47"-51"
	2X-Large	<b>2215.24</b>	50"-54"	52"-56"	51"-55"	53"-57"

**Women's Pendulous Abdomen with 15° Lordosis**

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2215.01PA</b>	28"-32"	29"-33"	26"-30"	35"-39"
	<b>2215.02PA</b>	30"-34"	31"-35"	29"-33"	37"-41"
	<b>2215.03PA</b>	32"-36"	34"-38"	32"-36"	38"-42"
	<b>2215.04PA</b>	36"-40"	37"-41"	37"-41"	41"-45"
	<b>2215.05PA</b>	41"-45"	42"-46"	43"-47"	45"-49"
	<b>2215.06PA</b>	43"-47"	46"-50"	45"-49"	50"-54"

**Women's Pendulous Abdomen with 0° Lordosis**

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2215.13PA</b>	28"-33"	30"-34"	30"-34"	35"-39"
	<b>2215.14PA</b>	31"-35"	32"-36"	32"-36"	36"-40"
	<b>2215.15PA</b>	33"-37"	35"-39"	35"-39"	39"-43"
	<b>2215.16PA</b>	37"-41"	38"-42"	38"-42"	41"-45"
	<b>2215.17PA</b>	42"-45"	43"-47"	43"-47"	45"-49"
	<b>2215.18PA</b>	44"-48"	48"-52"	49"-53"	51"-54"

**Men's Pendulous Abdomen with 15° Lordosis**

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2215.07PA</b>	34"-38"	34"-38"	32"-36"	35"-39"
	<b>2215.08PA</b>	36"-40"	36"-40"	35"-39"	38"-42"
	<b>2215.09PA</b>	39"-43"	38"-42"	38"-42"	41"-45"
	<b>2215.10PA</b>	42"-46"	42"-46"	42"-46"	44"-48"
	<b>2215.11PA</b>	45"-49"	47"-51"	47"-51"	48"-52"
	<b>2215.12PA</b>	51"-55"	51"-55"	51"-55"	53"-57"

**Men's Pendulous Abdomen with 0° Lordosis**

Quantity	Model No.	Axilla Circ.	Xyphoid Circ.	Waist Circ.	Troch Circ.
	<b>2215.19PA</b>	34"-38"	35"-39"	34"-38"	36"-40"
	<b>2215.20PA</b>	36"-40"	38"-42"	38"-42"	39"-43"
	<b>2215.21PA</b>	37"-41"	40"-44"	41"-45"	42"-46"
	<b>2215.22PA</b>	40"-44"	44"-48"	44"-48"	44"-48"
	<b>2215.23PA</b>	45"-49"	48"-52"	49"-53"	47"-51"
	<b>2215.24PA</b>	50"-54"	53"-56"	54"-58"	53"-57"



6333 North Orange Blossom Trail • Orlando FL 32810 • USA  
800-446-6770 | www.orthomerica.com

